

Associate

Dear Colleagues:

If you are interested in sending specimens to the laboratory at MERSI for O.C.T. processing for immunochemistry study, please do the following way (1 or 2):

- 1. Please be sure to provide us enough size of specimen (at least 3x5 mm²) with patient's brief clinical information.
 - 2. Immerse the specimen in the ZEUS fixative solution for no more than five days. Keep it at 4°C. Afterward, the specimen can be sent to us via express mail service in the ZEUS fixative solution. We will then embed it in Tissue-Tek O.C.T. compound for cryostat sectioning.
 - 3. You may get ZEUS fixative solution from: "ZEUS Scientific"

Cat. 0102

Tel: 1-800-286-2111 Option 2

Website: www.zeusscientific.com

- 4. Alternatively, the specimen can be snap-frozen with Liquid Nitrogen immediately after cutting the tissue. This frozen specimen can then be embedded in Tissue-Tek Compound (O.C.T.). It should always be kept at -80°C, and then it can be sent to us on **sufficient** dry ice by overnight mailing.
- 5. You may get (O.C.T.) from:
 "Richard-Allan Scientific"
 Name-Frozen section medium
 Cat. 6502
 Tel:800-522-7270
- Print Requisition Form Complete and mail with specimen

MERSI.COM | UVEITIS.ORG



Stephen D. Anesi, MD, FACS Partner and Co-President

Peter Y. Chang, MD, FACS

Partner and Co-President

Peter L. Lou, MD

Associate

7. Send to attn.: Dr. Tongzhen Zhao (to address above). Before sending out the specimen, please contact Dr. Zhao (email: tzhao@mersi.com) to notify her and confirm her availability. Our lab is closed on Saturdays and Sundays, so please refrain from sending it on Fridays.

If you have any questions, please contact Scott Evans at (781) 891-6377 or sevans@mersi.com.

Associate



Associate

Diagnostic Services Available.

1. **Routine Histopathology:** In addition to the standard Haemotoxylin and Eosin (H&E), Periodic acid-Schiff (PAS), and Giemsa staining, we use plastic embedding which is cut with gass knife to achieve very thin (0.1-0.2 micron) section of the ocular tissue.

2. Immunopathology

a.**Immunofluorescence:** This can be done for both conjunctival and scleral specimens. Staining is done for IgA, IgG, IgD, IgM, IgE, C3, C4, albumin (negative control), fibrinogen (positive control), and type IV collagen (positive control) to ensure high sensitivity and specificity in tissue diagnosis. In addition, immunofluorescent staining for anti-HSV and anti-VZV antibodies can be performed when a herpetic cause of conjunctivitis or scleritis is suspected.

b.**Immunoperoxidase:** This is only done on conjunctival specimen if the immunofluorescence result is equivocal.

Instructions for Specimen Submission: Please follow one of these procedures for tissue submission.

If you live in California, Florida, Maryland, New York, Pennsylvania, or Rhode Island, please call Scott Evans at 617-620-8033 before submitting specimens.

If your patient has the following insurance, please contact MERSI directly before sending specimen samples:

- Ambetter from CeltiCare
- BMC
- Celtic Care
- Commonwealth Care Alliance
- Empire NY
- Free Care HSN/Partial HSN

Massachusetts Eye Research and Surgery Institution

Partner and Co-President

Associate



Partner and Co-President

- FCHP
- GHI
- Golden Rule
- Kaiser Permanente
- Minuteman Health
- Optum Care
- Oxford Liberty
- International Insurance Plans
- Any form of Medicaid

For questions regarding insurance, please contact:

Ferida Tadic – MERSI Director of Operations: 781-891-6377

Ione Villafane – MERSI Supervisor of Billing and Surgical Scheduling: 781-891-6377



AssociateAssociate

MERSI IMMUNOPATHOLOGY LAB 1440 Main Street, Suite 201 Waltham, MA 02451

(Phone) 781-891-6377 (Fax) (781) 647-1430

TEST REQUISITION FORM

	PATIENT INFORMATION	PHYSICIAN INFORMATION
Patient Name:		Physician Name:
Patient ID#:		Institution:
DOB:	SEX:	Address:
Insurance Compar	ny Name:	
Insurance ID:		
Subscriber Name:	DOB:	
Self Pay: Circle On	e-Y N	
Credit Card Name	:	
Credit Card Numb	er:	
Credit Card Exp Da	ate:	
Lab Accession #:		
Collection Date:	Time:	Phone: FAX:
Specimen Type:		
Clinical History:		Physician Signature:
		Physician Name (printed):
Diagnosis:		Date Signed:

Associate

DIRECTOR: Dr. Stephen Lyle

Associate



	TESTS REQUESTED				
Tests		Comments			
	H&E				
	PAS				
	Giemsa				
	Immunofluorescence				
	Immunoperoxidase				
	Indirect Immunofluorescence				
SPECIMEN DELIVERY ADDRESS		CLINICAL LAB USE ONLY			
Send	samples OVERNIGHT on cold packs to:	Date received: Time:			
	MERSI IMMUNOPAYHOLOGY LAB	Specimen Type:			
ATT: Dr. Tongzhen Zhao		Report Date:			
1440 Main Street, Suite 201		Sample Condition:			
Waltham, MA 02451		Comments:			
Call lab (20 78.1-891-6377 prior to sending (ask for Dr. Zhao)See instructions for fixing tissue attached				
Ship samples only Monday through Thursday					

CLIA # 22D2017227